



# **New Age**

## **FARMA BIKE B-2**



Read the instructions carefully before use and keep this manual for future reference.



**New Age**

NEW AGE ITALIA S.R.L.  
Via De Brozzi, 3 – 48022 Lugo (RA)  
[www.newageitalia.it](http://www.newageitalia.it)

# ***FARMA BIKE B-2***

## ***PASSIVE MOTION DEVICE FOR REHABILITATION***

### **ADVANTAGES**

- Can help with the lack of exercise
- Train you wrist, leg and arm muscles.
- Helps strengthen your heart and lungs.
- Helps with blood circulation
- Slows down the aging process.

### **MAIN FEATURES**

Timer: 1-30 minutes

Speed 15 different speeds

### **CAUTION**

- If irregular heartbeat occurs during usage, stop immediately and contact your physician.
- Be extra careful of long hair and wool cloth for it might get caught in the machine.
- If the machine becomes damaged in anyway both appearance and mechanical do not use the machine. Unplug the machine and contact the store of purchase.
- During usage if strange odor occurs, or the machine begins to smoke or over heats, please unplug the machine.
- Please make sure that the wheel chair or chair that you are sitting on is properly secure before usage.
- Please do not use the machine barefooted. Wear proper shoes when using the machine.
- This machine is intended for one individual to use at one time. Do not attempt to use the machine with more two or more people.
- If the electrical is damaged in any way please do not use the machine in case it will cause an electrical fire.
- If there is water or water vapors in the electrical cord, please do not use the machine.
- Please do not use the machine immediately after a meal it might cause the user to become uncomfortable.
- Do not use the machine in high temperature or in a sauna, it might cause damage to the

machine.

- Please do not use a foreign object to strike the machine or drop the machine it will cause damage to the machine.
- Make sure to unplug the machine after usage.
- Please do not attempt to repair or alter the machine in any for it might cause the machine to malfunction.
- When usage please remember to set the timer, speed and settings. Do not over use the product for you might become uncomfortable.

## **PRODUCT CARE**

- Do not store the product in high temperature or heavy moisture.
- Always unplug the electrical cord.
- If the machine is dirty, please a basic cleaning solution to clean the machine. Do not use an acidic cleaning solution.
- Do not place heavy objects on the machine it might cause the machine to become deformed and cause malfunctions.
- Do not take apart the machine or try to modify it in any way.
- Over time if any damage occurs to the machine or the electrical wires, please take the machine to the store where it is purchased.

## **OPERATING THE MACHINE**

1. When using the product, please make use that the machine does not slide and is properly secured. It is suggested that the machine be use on an even ground.
2. Make sure that wheel chair or chair is properly secure in place. Make sure that it does not move or slide around for it might cause accidents.
3. If the place where you are using the machine is not on an even ground, please adjust the machine so that it is secure and safe to use,
4. The headrest can be adjusted to fit to a user's specific needs.
5. To adjust the machine to a user's needs, use the nut to extend or contract the machine.

## **HOW TO OPERATE THE AUTO-EXERCISE BIKE**


1. Please inset the power plug into the socket (panel will display the red light)
2. Please power on the handle and pedal button (panel will display the time)
3. A hand grip, double foot on the pedal board

4. Please use the up and down button to set the time (from 1 to 10 minutes)
5. Switch display speed button (the speed light will show up)
6. Please use the up and down button to adjust the speed you want (from 1 to 15 grades)

**Caution:** The user must put their both feet on the pedal when set up the speed and the speed must to match the user's situation. Excessive speed may cause harm.

- The distance and calories will be showed on the panel.
- The distance is from 1-10 km
- When the setting time is up, the machine will automatically turn off. All the data will turn to zero.
- Turn off the machine before leaving.
- Time, distance and calories are for reference. The actual data will be a little bit different.

## SPECIFICATIONS



1. Material: Steel and ABS
2. Grip and pedal: Polyurethane
3. Dimension: approx. 70 x 58 x 81-101 (high)cm
4. Power: 220V/50W
5. Timer: 1-30 minutes
6. Rotation frequency: approx. 18 to 80 per minute
7. Wire length: approx. 2.5 m
8. Weight: approx. 15 kg
9. Automatic Stop Protection System
10. Close timer feature

For assistance, please contact:

**NEW AGE ITALIA S.R.L.**

Via De Brozzi, 3 – 48022 Lugo (RA)

[www.newageitalia.it](http://www.newageitalia.it) [info@newageitalia.it](mailto:info@newageitalia.it)

TEL +39 0545 32019

FAX +39 0545 369028