



LASERTRON CARD

LASERTRON Card is the new professional instrument for lasertherapy, characterized by the numerous functions and possibility of use. With the head in endowment, the instrument can disburse the laser light to effect numerous therapeutic applications. LASERTRON Card contains pre-setted programs with the most common parameters for the applications and a free space in memory to plan new personalized programs. LASERTRON Card can be fed by electric current or with the rechargeable inside battery. It is endowed with a comfortable handbag for the transport and the home therapy.



PROGRAMS/APPLICATIONS

- Inflammations
- Tensond pain
- Bursitises
- Bruises
- Muscular tears
- Distorsions
- Articular pains
- Arthritises and arthrosis
- Cellulitis
- Scars
- Problems of the skin
- Acute and chronic pains



TECHNICAL CARACHTERISTICS

Power suply :battery, external feeder

Output channels: 2 independent

Programs: 20 pre-setted + 20 free

Laser source 1: pulsated diode Ga Á the, power 38mW (60W peak), wavelength 903nm

Laser source 2: pulsated diode Ga As IR, power 2mW, wavelength 903nm

Frequency : 10/10.000Hz

Radiate issue: continuous, pulsated

Dimensions 18*17*10cm, 1 kg only device, 2,5 kg with and weight accessories

BASIS ENDOWMENT

- Suitcase
- 2 laser probes
- Rechargeable battery
- External feeder
- Manual of use